



<b>Fecal Tests</b>	<b>Panel includes:</b>	<b>Tube/specimen requirements</b>	<b>Turn around time</b>
<b>Fecal <i>Giardia</i></b>	Fecal <i>Giardia</i> testing by ELISA methods	1-3 grams of fresh feces submitted in a sterile container (fecalyzer). <b>If fecal float is required submit two separate samples. Samples will not be saved or split.</b>	Run weekly on Tuesdays and Fridays in the morning. Results are available after 2pm.
<b>Fecal Ova and Parasite Float</b>	Fecal flotation by zinc sulfate centrifugation	1-3 grams of fresh feces submitted in a sterile container (fecalyzer). <b>If fecal <i>Giardia</i> is requested submit two separate samples. Samples will not be saved or split.</b>	1 day